Calcitonin II

¹²⁵I RIA Kit

For the quantitative determination
of human calcitonin in serum

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1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE.

This kit contains instructions and materials for the quantitative determination of human calcitonin in serum by radioimmunoassay (RIA).

2. SUMMARY AND EXPLANATION

Calcitonin is a peptide hormone with a molecular weight of 3,418 and contains thirty two amino acids:

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17
H-Cys-Gly- Asn-Leu- Ser- Thr- Cys-Met- Leu- Gly- Thr- Tyr- Thr- Gln- Asp-Phe- Asn-
18 19 20 21 22 23 24 25 26 27 28 29 30 31 32
Lys- Phe- His- Thr- Phe-Pro- Gln- Thr- Ala- Ile- Gly- Val- Gly- Ala- Pro- NH₂

In humans, this hypocalcemic hormone is secreted by thyroidal parafollicular cells of neuroectodermal origin, probably in response to hypercalcemia.^{10,14} Calcitonin's hypocalcemic actions are mediated via effects on bone and kidney.⁶

DiaSorin has developed and thoroughly tested a sensitive calcitonin radioimmunoassay (RIA) that can detect calcitonin (CT) concentrations as low as 15 pg/mL. The primary use for calcitonin RIA is serological detection of malignant growths. Medullary thyroid cancer (MTC), which comprises 4-8% of thyroid tumors, secretes calcitonin in excess. Very high levels of circulating calcitonin, detected by RIA, are useful in evaluating this type of tumor.^{5,7,11,21} Some of these MTCs are hereditary, and it is important to screen **all** members of a family for calcitonin excess when a case is discovered. Other uses of the calcitonin RIA include investigation of ectopic calcitonin production in cancer patients, particularly those with breast and lung tumors.^{4,9,13,18,24} Very active investigation is taking place to discover the usefulness of calcitonin measurements in the evaluation of nonthyroidal type malignancies.^{8,16,19}

3. PRINCIPLES OF THE ASSAY

The DiaSorin Calcitonin II RIA is a disequilibrium procedure using delayed tracer addition to increase sensitivity. The antibody was produced in a goat against pure synthetic human calcitonin. In this RIA, sample and first antibody are combined and incubated for 16-24 hours at 2-8°C. Tracer is then added, followed by a second incubation for 16-24 hours at 2-8°C. A pre-precipitated second antibody complex is added to separate the bound from free tracer. The assay can then be centrifuged and decanted after 15-20 minutes incubation at 20-25°C.

4. REAGENTS PROVIDED IN THE KIT

Calcitonin Calibrator 0	1 vial/10 mL	2 vial/10 mL
Calcitonin Calibrator	2 vials/1.0 mL	4 vials/1.0 mL
Calcitonin Antiserum	1 vial/14 mL	2 vials/14 mL
¹²⁵ I Calcitonin	1 vial/7 mL	2 vials/7 mL
Calcitonin Precipitating Complex	1 vial/35 mL	2 vials/ 35 mL
Calcitonin Controls	2 vials/1.0 mL	4 vials/ 1.0 mL
Number of tests	65	130

STORAGE: Upon receipt, and prior to reconstitution, store all reagents at 2-8°C. After reconstitution store all reagents at -15° or lower until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

4.1 Calcitonin Calibrator 0: lyophilized reagent

BSA-borate buffer with sodium azide (0.25%) added. Reconstitute the vial with 10 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.2 Calcitonin II Calibrator: lyophilized reagent

Human synthetic calcitonin, at a nominal concentration of 1,000 pg/mL, is diluted in BSA-borate buffer with sodium azide (0.2%) and other stabilizers added. Exact concentration values are assigned with each lot. Reconstitute the vial with 1.0 mL of purified water, mix and allow it to stand for 15-20 minutes on crushed ice until the contents are completely dissolved; mix thoroughly before using. In order to obtain the entire calibrator curve, make serial dilutions by adding 500 µL of calibrator to 500 µL of calibrator 0. For example, if the calibrator concentration is 1,000 pg/mL, calibrators of 500, 250, 125, 62.5 and 31.25 pg/mL will be made as follows:

Add 500 µL of 1,000 pg/mL calibrator to 500 µL of calibrator 0 and mix to give 500 pg/mL.

Add 500 µL of 500 pg/mL calibrator to 500 µL of calibrator 0 and mix to give 250 pg/mL.

Add 500 µL of 250 pg/mL calibrator to 500 µL of calibrator 0 and mix to give 125 pg/mL.

Add 500 µL of 125 pg/mL calibrator to 500 µL of calibrator 0 and mix to give 62.5 pg/mL.

Add 500 µL of 62.5 pg/mL calibrator to 500 µL of calibrator 0 and mix to give 31.25 pg/mL.

CAUTION: Keep tubes on crushed ice while the serial dilutions are being made and prior to use in the assay. The DiaSorin calcitonin calibrator has been analyzed against the World Health Organization calcitonin standard 70/234. The calibrator demonstrates commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

DISCARD ANY EXCESS; DO NOT FREEZE.

4.3 Calcitonin Antiserum: lyophilized reagent

Goat anti-calcitonin serum is diluted in BSA-borate buffer with sodium azide (0.1%) and blue dye added. Reconstitute the vial with 14 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.4 ¹²⁵I Calcitonin: lyophilized reagent

Synthetic human calcitonin is labeled with iodine-125 and diluted in BSA-borate-EDTA buffer with sodium azide (0.4%) and red dye added. Reconstitute the vial with 7 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix gently before using.

4.5 Calcitonin Precipitating Complex: lyophilized reagent

Normal goat serum, pre-precipitated with donkey anti-goat serum and polyethylene glycol (PEG), is diluted in BSA-borate buffer with sodium azide (0.1%) added. Reconstitute the vial with 35 mL of purified water. Mix thoroughly until the suspension appears homogenous and then allow it to stand for a minimum of thirty minutes at room temperature with occasional mixing.

4.6 Calcitonin Controls (Level 1 and 2): lyophilized reagent

Human serum is spiked, if necessary, with the appropriate amount of synthetic human calcitonin to obtain a concentration within a specified range. Sodium azide (0.1%) and other stabilizers are added. Reconstitute the vial with 1.0 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly and treat as an unknown sample. Ranges are printed on the control vials.

CAUTION: DISCARD ANY EXCESS; DO NOT FREEZE.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 1.5 μ Ci (55.5 kBq) kit No. 25065 or 3 μ Ci (111 kBq) kit No. 25130 of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the

package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERORATION OF KIT REAGENTS

- 6.1** The presence of abnormal particulate matter in any of the reagents.
- 6.2** A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3** A decrease in maximum binding.
- 6.4** A high nonspecific binding.

7. COLLECTION AND PREPARATION OF THE SPECIMEN

One hundred microliters, in duplicate, of serum are required for the assay. Collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. Centrifuge for fifteen minutes using 760 x g* to obtain hemolysis-free serum. No additives or preservatives are required to maintain integrity of the sample. All plastics, glassware or other materials coming into contact with the specimen should be entirely free of any contamination. A fasting specimen is recommended but not required. Calcitonin is labile at room temperature. If the serum is not assayed immediately, it should be stored at -15°C or lower. Specimens should not be repeatedly frozen and thawed.

8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1** Disposable borosilicate glass tubes, 12 x 75 mm. Plastic tubes **cannot** be used.
- 8.2** Temperature controlled centrifuge to accommodate 12 x 75 mm tubes (20-25°C).
- 8.3** Gamma scintillation counter capable of counting iodine-125.
- 8.4** Vortex.
- 8.5** Pipetting devices:
 - a.** Micropipettors calibrated to deliver 100 µL and 200 µL.
 - b.** Repeating dispensers calibrated to deliver 200 µL and 500 µL.
- 8.6** Purified water

9. ASSAY PROCEDURE

- 9.1** Thaw unknown samples and place on crushed ice.
- 9.2** Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 20-25°C. Mix gently and then place on ice before using.
- 9.3** Set up labeled 12 x 75 mm disposable borosilicate glass tubes in duplicate according to the Scheme of the Assay, on the back page.
- 9.4 Place the rack of tubes on crushed ice.**
- 9.5** Add reagents as follows:
 - a. Total count tubes**
Set aside until step 7
 - b. Nonspecific binding (NSB)**
100 µL of calibrator 0
 - c. Calibrator 0**
100 µL of calibrator 0
200 µL of calcitonin antiserum (blue)
 - d. Calcitonin calibrators**
100 µL of calcitonin calibrator
200 µL of calcitonin antiserum (blue)

e. Quality control and unknown samples

100 µL of serum

200 µL of calcitonin antiserum (blue)

- 9.6** Vortex the tubes gently without foaming and incubate for 16-24 hours at 2-8°C.
- 9.7** Add 100 µL of ¹²⁵I calcitonin (red) to all tubes.
- 9.8** Vortex the tubes gently without foaming and incubate for 16-24 hours at 2-8°C.
- 9.9** Vigorously mix the precipitating complex; add 500 µL to all tubes except the total count tubes.
- 9.10** Vortex the tubes gently without foaming and incubate for 15-25 minutes at 20-25°C.
- 9.11** Centrifuge using 760 x g* for twenty minutes at 20-25°C.
- 9.12** Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a maximum of two minutes. Blot the tubes with absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
- 9.13** Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy. (See Limitations of the Procedure section.)

10. PROCEDURAL COMMENTS

- 10.1** Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
- 10.2** Assay tubes must be set up on crushed ice to avoid spuriously high values.
- 10.3** Some manufacturers' disposable borosilicate glass tubes yield elevated nonspecific bindings. **Do not** use plastic tubes in this assay.
- 10.4** If you choose to aspirate supernatant from precipitate, be careful not to disturb the precipitate.
- 10.5** To completely monitor the consistent performance of an RIA there are additional factors which may be checked. DiaSorin suggests a check of the following parameters to assure consistent kit performance.
- a. Total Counts**
 - b. Maximum Binding**
Average counts per minute (CPM) of calibrator 0 Tubes / Average CPM of Total Count Tubes.
 - c. Nonspecific Binding**
Average CPM of NSB Tubes / Average CPM of Total Count Tubes.
 - d. Slope of Calibrator Curve**
For example, monitor the 80, 50 and 20% suppression points of the calibrator line.

11. QUALITY CONTROL

Each laboratory should include at least two control sera in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten assays. An acceptable range of values may then be obtained for these controls using ± 2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the quality control sera included in this kit.

$$* \quad g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

12. CALCULATION OF RESULTS

There are many methods in existence for calculating results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either a linear or logarithmic scale. Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin Quality Control Laboratory is % B/B₀ versus log concentration.

12.1 Calculate the average CPM for each calibrator, control and unknown sample.

12.2 Subtract the average CPM of the NSB tubes from all counts.

12.3 Divide the corrected CPM of each calibrator, control or unknown sample by the corrected CPM of the calibrator 0.

$$B/B_0 \% = \frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of 0 Calibrator} - \text{CPM of NSB}} \times 100$$

12.4 Using 2 cycle semi-log or log-logit graph paper, plot percent B/B₀ for the calcitonin calibrators (vertical axis) versus the concentration (horizontal axis).

12.5 Draw a best-fit line through the points.

12.6 Interpolate the levels of calcitonin in the unknown samples from the plot.

12.7 If any unknown sample reads greater than the highest calibrator, it should be diluted appropriately with calibrator 0 and assayed again.

12.8 If an unknown sample has been diluted, correct for the appropriate dilution factor.

12.9 Calculate maximum binding by dividing CPM of calibrator 0 by the average total counts obtained in the total count tubes.

TABLE II DiaSorin Calcitonin II RIA Sample Data					
Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent Bound (B/T)	Percent Conc. (B/B ₀)(pg/mL)
Total Count	17,943 18,128	18,036			
NSB	878 824	851		4.7	
Calibrator 0	8,055 8,049	8,052	7,201	39.9	100.0
Calibrators (pg/mL)					
A (35.0)	7,344 7,211	7,278	6,427		89.3
B (70.0)	6,969 6,649	6,809	5,958		82.8
C (140)	5,390 5,476	5,433	4,582		63.7
D (280)	3,628 3,725	3,676	2,826		39.3
E (560)	2,301 2,286	2,293	1,442		20.0
F (1,120)	1,519 1,546	1,532	682		9.5
Unknown Samples					
1	6,087 5,940	6,014	5,162		71.7 115
2	4,796 5,012	4,904	4,053		56.3 163

Typical sample data and a calibrator curve are shown in TABLE II and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.


FIGURE 1

13. LIMITATIONS OF THE PROCEDURE

- 13.1** If the initial concentration of any unknown sample is greater than the highest calibrator, dilute with calibrator 0 only.
- 13.2** Counting times should be sufficient to prevent statistical error (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).
- 13.3** The DiaSorin QC lab uses a smoothed spline curve fit.

14. EXPECTED VALUES

Normal Range

Each laboratory should establish its own normal range. Normal values observed in the DiaSorin Laboratory were 47 (± 24 pg/mL). This would lead to a normal range (2 std. dev.) of nondetectable (ND)-95. Smoking may cause slightly elevated serum calcitonin levels.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Precision

Average Within Assay Variation

	Mean Value (pg/mL)	S.D	%C.V.
Low	67.3	10.1	14.8
Medium	278.5	18.9	6.7
High	530	30.7	5.7

Between Assay Variation

	Mean Value (pg/mL)	S.D	%C.V.
Low	92.76	13.42	14.47
Medium	142.55	16.55	11.61
High	528.79	36.7	6.95

15.2 TRUENESS: THE TRUENESS OF THE ASSAY HAS BEEN VERIFIED BY THE LINEARITY TEST AND THE RECOVERY TEST.

Linearity (Parallelism)

Serial Dilution Study of four Unknown Samples

Sample Number	Undiluted (pg/mL)	1/2	1/4	1/8
1	534	624.5	663.7	696.1
2	277.2	291.8	346.5	289.4
3	975.6	1216.9	1396.8	1487.6
4	819.3	1047.9	1144.2	1081

Recovery

Five unknown samples have been diluted 1:1 with calibrators 0, 48, 96, 193, 385, 770 and 1,540 pg/mL. Each dilution was then assayed and an overall mean recovery % was calculated for every unknown sample used.

Sample Number	Overall Mean Recovery (%)
1	105
2	83
3	77
4	85
5	85

15.3 Analytical Sensitivity

When defined as the apparent concentration at three standard deviations from the counts at maximum binding, the minimum detectable amount is 1.5 pg/tube (15 pg/mL).

15.4 Analytical Specificity

Comparison of the cross-reactivity of calcitonin antibody was made with the following peptides:

Peptide	% Cross-reactivity
Parathyroid Hormone	< 0.01
Insulin	< 0.01
Thyroid Stimulating Hormone	< 0.01
Adrenocorticotrophic Hormone	< 0.01
Prolactin	< 0.01

Salmon calcitonin does not cross-react with the antisera. Rat calcitonin does cross-react and the assay can be used for studies on rats.

15.5 Interference

Interference studies were performed using NCCLS EP7-A as a guideline, to assess the effects of common endogenous interferents. Values obtained from demonstrated interference due to the addition of Triglyceride and hemoglobin. The addition of cholesterol and Bilirubin did not effect results.

SEE LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow reagents to reach above 25°C. Mix gently then place on crushed ice before using.
2. Identify tubes in duplicate. Place the rack of tubes on crushed ice.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Calibrator 0	-	100 µl	100 µL	-
Calibrators (1-6)	-	-	100 µL	-
Extracted Controls	-	-	-	100 µL
Extracted Unknown	-	-	-	100 µL
Samples	-	-	200 µL	200 µL
Antiserum				

4. Cover the tubes with Parafilm and vortex gently.
5. Incubate for 16–24 hours at 2–8°C.
6. Dispense 100 µL of Tracer into all wells.
7. Cover and vortex gently.
8. Incubate for 16–24 hours at 2–8°C.
9. Dispense 500 µL of precipitating complex into all wells except the Total Count tubes.
10. Incubate for 15–25 minutes at 20–25°C.
11. Centrifuge using 760 x g* for 20 minutes.
12. Decant the supernatants.
13. Count each tube in a gamma counter for 60 seconds or longer.

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
REFERENCES/BIBLOGRAPHIE/LITERATUR/BIBLIOGRAFÍA/BIBLIOGRAFI

A

1. Abe, Kaoru, I. Adachi, S. Miyakawa, M. Tanaka, K. Yamaguchi, N. Tanaka, T. Kameya and Y. Shimosato, "Production of Calcitonin, Adrenocorticotrophic Hormone and β -Melanocyte-stimulating Hormone in Tumors Derived from Amine Precursor Uptake and Decarboxylation Cells," **Cancer Research**, 37:4190, (1977).
2. Baylin, S.B., W.R. Weisberger, J.C. Eggleston, G. Mendelsohn, M.A. Beaven, M.D. Abeloff and D.S. Ettinger, "Variable Content of Histaminase, L-Dopa Decarboxylase and Calcitonin in Small-Cell Carcinoma of the Lung," **The New England Journal of Medicine**, 299(3):105, (1978).
3. Bertagna, X.Y., W.E. Nicholson, O.S. Pettengill, D.G. Sorenson, C.D. Mount and D.M. Orth, "Ectopic Production of High Molecular Weight Calcitonin and Corticotropin by Human Small Cell Carcinoma Cells in Tissue Culture: Evidence for Separate Precursors," **The Journal of Clinical Endocrinology and Metabolism**, 47(6):1390, (1978).
4. Coombes, R.D., G.C. Easty, S.I. Detre, C.J. Hillyard, U. Stevens, S.I. Girgis, L.S. Galante, L. Heywood, I. MacIntyre and A.M. Neville, "Secretion of Immunoreactive Calcitonin by Human Breast Carcinomas," **British Medical Journal**, 4:197, (1975).
5. Corwin, T.R., "Medullary Carcinoma of the Thyroid," **Surgery, Gynecology and Obstetrics**, 138:453, (1974).
6. Deftos, L.B., B. Roos and J. Parthemore, "Calcium and Skeletal Metabolism," **Western Journal of Medicine**, 123:447, (1975).
7. Goltzman, D., J.T. Potts, Jr., E.C. Ridgeway and F. Maloof, "Calcitonin as a Tumor Marker – Use of the Radioimmunoassay for Calcitonin in the Post-operative Evaluation of Patients with Medullary Thyroid Carcinoma," **The New England Journal of Medicine**, 290(19):1035, (1974).
8. Gropp, C., K. Havemann and A. Scheuer, "Ectopic Hormones in Lung Cancer Patients at Diagnosis and During Therapy," **Cancer**, 46:347, (1980).
9. Hansen, M., H.H. Hansen and N. Tryding, "Small Cell Carcinoma of the Lung: Serum Calcitonin and Serum Histaminase (Diamine Oxidase) at Basal Levels and Stimulated by Pentagastrin," **Acta Medica Scandinavica**, 204:257, (1978).
10. Heath, III, J. and G.W. Sizemore, "Plasma Calcitonin in Normal Man – Differences Between Men and Women," **The Journal of Clinical Investigation**, 50:1135, (1977).
11. Hennessy, J.F., T.K. Gray, C.W. Cooper and D.A. Ontjes, "Stimulation of Thyro-calcitonin Secretion by Pentagastrin and Calcium in Two Patients with Medullary Carcinoma of the Thyroid," **The Journal of Clinical Endocrinology**, 7:45, (1977).
12. Hill, Jr., C.S., "Medullary Carcinoma of the Thyroid," **American Federation Proceedings**, 99, (1973).
13. Hillyard, C.J., D.G. Oscier, R. Foa, D. Catovsky and J.M. Goldman, "Immunoreactive Calcitonin in Leukaemia," **British Medical Journal**, 2:1392, (1979).
14. Parthemore, J.G. and L.J. Deftos, "Calcitonin Secretion in Normal Human Subjects," **Journal of Clinical Endocrinology and Metabolism**, 47(1):184, (1978).
15. Rees, L.H., "Ectopic Hormone Production by Cancer Cells," **Lab-Lore** 8, (1978).
16. Richardson, R.L., F.A. Greco, R.K. Oldham and G.W. Liddle, "Tumor Products and Potential Markers in Small Cell Lung Cancer," **Seminars in Oncology**, 5:253, (1978).
17. Roof, B.S., R. Weinstein, I., Vuji, N.M. Burdash and M. Dopson, "Immunoheterogeneity of the Calcitonins of Hypercalcemia, Breast and Lung Cancer and Medullary Carcinoma of Thyroid," **Abstract 797**, Endocrine Society Meeting No. 60, Miami, Florida, (1978).


18. Roos, B.A., "Calcitonin as Tumor Marker," **RIA for Physicians** 5, (1980).
19. Roos, B.A., A.W. Lindall, S.B. Baylin, J. O'Neil, A.L. Frelinger, R.S. Birnbaum and P.W. Lambert, "Plasma Immunoreactive Calcitonin in Lung Cancer," **Endocrine Research Communications**, 6(2):169, (1979).
20. Roos, B.A., A.W. Lindall, S.B. Baylin, J. O'Neil, A.L. Frelinger, R.S. Birnbaum and P.W. Lambert, "Plasma Immunoreactive Calcitonin in Lung Cancer," **Journal of Clinical Endocrinology and Metabolism**, 50:659, (1980).
21. Samaan, J.A., M. Ibanez and C.S. Hill, Jr., "Medullary Carcinoma of the Thyroid and Astrocytoma," **Annals of Internal Medicine**, 86(5):585, (1977).
22. Schwartz, K.E., A.R. Wolfson, B. Forster and W.D. Odell, "Calcitonin in Nonthyroidal Cancer," **Journal of Clinical Endocrinology and Metabolism**, 49:438, (1979).
23. Schwartz, K.E., A.R. Wolfson and W.D. Odell, "Plasma Calcitonin – Use as a Tumor Marker," **Abstract 572**, Endocrine Society Meeting No. 60, Miami, Florida, (1978).
24. Silva, O.L. and K.L. Becker, "High Plasma Calcitonin Levels in Breast Cancer," **British Medical Journal**, 21:460, (1976).
25. Silva, O.L., L.E. Broder, J.L. Doppman, R.H. Snider, C.F. Moore, M.H. Cohen and K.L. Becker, "Calcitonin as a Marker for Bronchogenic Cancer," **Cancer**, 44:680, (1979).
26. Silva, O.L., K.L. Becker, A. Primack, J.L. Doppman and R.H. Snider, "Hypercalcitemia in Bronchogenic Cancer—Evidence for Thyroid Origin of the Hormone," **Journal of the American Medical Association**, 234(2):183, (1975).
27. Silva, O.L., K.L. Becker, A. Primack, J.L. Doppman and R.H. Snider, "Increased Serum Calcitonin Levels in Bronchogenic Cancer," **Chest**, 69(4):495, (1976).
28. Tashjian, A.H. Jr., B.G. Howland, K.E.W. Melvin and C.S. Hill, Jr., "Immunoassay of Human Calcitonin—Clinical Measurement, Relation to Serum Calcium and Studies in Patients with Medullary Carcinoma," **The New England Journal of Medicine**, 283 (17):890, (1978).

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Antiserum	Antisérum	Antiserum	Antisero	Antisiero
	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reactivo precipitante	Reagente precipitante
	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Trazador: antígeno etiquetado con ¹²⁵ I	Tracciatore: antigene etichettato con ¹²⁵ I
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Radioactive	Radioactif	Radioaktiv	Radiactivo	Radioattivo
	Harmful	Nocif	Gesundheits-schädlich	Nocivo	Nocivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

Continued from previous page - Symbols used with devices

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Phone: 651-439-9710
FAX: 651-779-7847

For Customer Service in the
U.S. and Canada Call Toll Free:
800-328-1482

In the United Kingdom Call:
44 118 9364200
FAX: 44 118 9792061

10399

27305
9/03

PRINTED IN U.S.A.

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DiaSorin S.A.
DiaSorin Inc.
Stillwater, MN 55082-0285 U.S.A.